



Metal Syringe

Weiyou® Veterinary metal syringe is made of high quality copper and stainless steel. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.



Product Description

Weiyou® Veterinary metal syringe is made of high quality copper and stainless steel. Weiyou® Veterinary metal syringe is mainly used for Cattle, Horse Sheep, Goat and other animals. OEM customization is acceptable. Vaccination to prevent disease for all kinds of animals. Weiyou® Veterinary metal syringe is continuous injection and accurate scale. Free samples, and small quantity can be supplied.

Specification

Item	Veterinary metal syringe
Usage	Vaccination
Material	Copper, stainless steel
MOQ	100pcs
Specification	5ml, 10ml, 20ml, 30ml, 50ml
Use For	Cattle, Horse Sheep, Goat, etc
Advantage	Continuous injection, accurate scale
Type	Vaccination to prevent disease
Packing	Carton

Features

Luer-lock adaptor

Brass raw material with chrome plated



www.nbweiyou.com

Glass tube easy for cleaning
Adjustable rubber piston rod









